

## ***elma lab clean A25 (ELC A25)***

STRONGLY ALKALINE FOAM-INHIBITED CLEANING CONCENTRATE FOR  
DIP & SPLASH CLEANING OF INSTRUMENTS IN LAB'S & WORKSHOPS

### **Description**

Suitable for labware and medical instruments from metal, glass, ceramic and plastic materials. Not for aluminium, magnesium, alloys of light-metals and alkali-sensitive glassware.

Removes markings, labels, greases from ground-glass joints; resinified & tarry residues; lime soap deposits and prevents from its redeposition; blood, saliva; protein-, bone- & tissue-residues; fats/greases and oils, polishing & grinding residues; fingerprints and dust. For strong/thick tarry residues the alkalinity of the cleaning solution may be amplified by the addition of potassium hydroxide or sodium hydroxide.

### **Application and dosage**

Ultrasonic bath: 1 - 2 vol% in water, at 50 - 75°C. Splash cleaning: 0.5 - 1 vol%, well above 55°C only. Tap or deionised water can be used preparing the cleaning solution.

Removal of *still fresh* blood, saliva or protein residues below 42°C recommended. Milkyness of heated cleaning solution does not reduce the cleaning performance. Rinse the parts after the cleaning and then dry.

### **Safety recommendations**

*elma lab clean A25 (ELC A25)* is classified as hazardous according to the regulation (EC) No 1272/2008 [GHS] (skin corrosive and serious eye damage, corrosive to metals [aluminium and light metals]). Observe also with respect to this the hints indicated in the Safety Data Sheet and always handle chemicals with care.

### **Physical-chemical Characterisation**

Density: ~1,14 g/ccm, pH (1 % aqueous solution): 12. Ingredients according to Annex VII A, EC-Regulation 648/2004 (detergents): 5-15% amphoteric surfactants, 5-15% non-ionic surfactants, <5% phosphates, <5% polycarboxylates.

### **Disposal**

The surfactants in our product meet the criteria for biodegradation as laid down in Annex III of the Regulation (EC) No 648/2004 on detergents.

The cleaning bath can be fed into the public sewage system after neutralisation; observe the local pH limit values and make sure that the contamination contained complies with the local sewage regulations. Neutralisation: acetic or citric acid in an ultrasonic bath – do not use hydrochloric or sulphuric acid there.

European waste code: 20 01 29\*, „detergents containing dangerous substances“.

### **Volumes, storage and transport**

Available volumes: 1 litre PE-bottle; 2.5 litre, 10 litre and 25 litre HDPE-can.

Store in closed original container at a temperature between +5°C and +30°C, protected from heat and direct solar radiation. After subcooling shake before use. Do not store with acids.

Shelf life: 3 years from date of production (see stamp on label).

Classification for all means of transport: class 8, UN 1814.

### **Elma Schmidbauer GmbH**

Gottlieb-Daimler-Str. 17, D-78224 Singen

E-Mail: [info@elma-ultrasonic.com](mailto:info@elma-ultrasonic.com)

15.09.2016 / (GB) Version 1.1

Tel.: +49 7731 882-0

Fax: +49 7731 882-266

[www.elma-ultrasonic.com](http://www.elma-ultrasonic.com)

